

August 2022

Dear valued customer,

Thank you for your inquiry regarding the expected service life of Stryker's Powered System.

The Power-PRO™ 2 powered ambulance cot, Power-PRO XT powered ambulance cot, Power-LOAD® powered fastener and lithium-ion battery charger all have a seven-year expected service life under normal use and with appropriate periodic maintenance.

The lithium-ion battery has a two-year expected service life following the same guidelines. End of battery service life will be indicated by a constant yellow LED to the far right of the battery charge indicator.

Stryker's products have been designed for their respective life expectancies under normal use conditions and with appropriate periodic maintenance as described in the maintenance manual for each device. Products in the field may perform beyond the expected useful life for several reasons, including if the device is not used as much as anticipated.

Feel free to contact us with any questions.

Sincerely,



John Guyeskey
Senior Marketing Manager

Copyright © 2022 Stryker.
M0000010681 REV AA

Stryker or its affiliated entities own, use, or have applied for the following trademarks or services marks: Power-PRO, Power-LOAD, Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.